



Product description

Order number	70.1183.101
Product description	Pipette tip, working volume: 5 ml, transparent, filling level rings, suitable for Gilson, 50 piece(s)/box

Product characteristics

Graduation	yes
Suitable for	Gilson
Volume of work	5 ml
Filter	no
Autoclavable (max. °C)	121 °C

Size

Max. volume	5 ml
Length of product	123 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 500
Type of inner packaging box
Piece(s) / inner box 50
Piece(s) / case 500
Piece(s) / pallet 16500
Case volume 0,0282 cbm
Weight of product 0,0065 kg
Weight of case 3,7 kg
EAN of inner box 4038917126465
EAN case 4038917080460

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

