



### **Product description**

Order number	70.1183.102
Product description	Pipette tip, working volume: 5 ml, transparent, filling level rings, suitable for Eppendorf, 50 piece(s)/box

### **Product characteristics**

Graduation	yes
Suitable for	Eppendorf
Volume of work	5 ml
Filter	no
Autoclavable (max. °C)	121 °C

### **Size**

Max. volume	5 ml
Length of product	123 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

**Purity & certification**

Product category	no medical device   IVD
Batched	yes

**Packaging**

Minimum order qty.	500
Type of smallest subpackaging	box
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	16500
Depth of case	165 mm
Width of case	311 mm
Height of case	610 mm
Case volume	0.0313 cbm
Weight of product	0.0065 kg
Weight of case	3.7 kg
EAN of inner box	4038917126472
EAN case	4038917080477

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

