



### **Product description**

Order number	70.1183.102
Product description	Pipette tip, working volume: 5 ml, transparent, filling level rings, suitable for Eppendorf, 50 piece(s)/box

### **Product characteristics**

Graduation	yes
Suitable for	Eppendorf
Volume of work	5 ml
Filter	no
Autoclavable (max. °C)	121 °C

### **Size**

Max. volume	5 ml
Length of product	123 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Material & colours**

Product material Polypropylene (PP)  
Colour of product transparent

### **Purity & certification**

Product category no medical device | IVD  
Batched yes

### **Packaging**

Minimum order qty. 500  
Type of inner packaging box  
Piece(s) / inner box 50  
Piece(s) / case 500  
Piece(s) / pallet 16500  
Case volume 0,0282 cbm  
Weight of product 0,0065 kg  
Weight of case 3,7 kg  
EAN of inner box 4038917126472  
EAN case 4038917080477

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

