



### Product description

|                     |   |
|---------------------|---|
| Order number        | 70.1186.100   |
| Product description | Pipette tip, working volume: 1,250 µl, transparent, filling level rings, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson and Brand and products of identical design, 96 piece(s)/box |

### Product characteristics

|                        |  |
|------------------------|--|
| Graduation             | yes  |
| Suitable for           | SARSTEDT Sarpette® M, Eppendorf, Gilson and Brand and products of identical design |
| Volume of work         | 1.250 µl   |
| Filter                 | no   |
| Coding colour          | blue   |
| Autoclavable (max. °C) | 121 °C   |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Size**

|                   |         |
|-------------------|---------|
| Max. volume       | 1250 µl |
| Length of product | 103 mm  |

**Material & colours**

|                   |                    |
|-------------------|--------------------|
| Product material  | Polypropylene (PP) |
| Colour of product | transparent        |

**Purity & certification**

|                  |                         |
|------------------|-------------------------|
| Product category | no medical device   IVD |
| Batched          | yes                     |

**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 3840          |
| Type of smallest subpackaging | box           |
| Piece(s) / inner box          | 96            |
| Piece (s) / outer case        | 3840          |
| Piece(s) / pallet             | 46080         |
| Depth of case                 | 264 mm        |
| Width of case                 | 475 mm        |
| Height of case                | 500 mm        |
| Case volume                   | 0.0627 cbm    |
| Weight of product             | 0.0022 kg     |
| Weight of case                | 9.05 kg       |
| EAN of inner box              | 4038917126489 |
| EAN case                      | 4038917080347 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

