



Product description

Order number	70.1186.100
Product description	Pipette tip, working volume: 1,250 µl, transparent, filling level rings, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson and Brand and products of identical design, 96 piece(s)/box

Product characteristics

Graduation	yes
Suitable for	SARSTEDT Sarpette® M, Eppendorf, Gilson and Brand and products of identical design
Volume of work	1.250 µl
Filter	no
Coding colour	blue
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Size

Max. volume 1250 µl
Length of product 103 mm

Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Product category not regulated
CE certified no CE mark
Batched yes

Packaging

Minimum order qty. 3840
Type of smallest box
subpackaging
Piece(s) / inner box 96
Piece (s) / outer case 3840
Piece(s) / pallet 46080
Depth of box 120 mm
Width of box 85 mm
Height of box 114 mm
Depth of case 495 mm
Width of case 267 mm
Height of case 466 mm
Case volume 0.0616 cbm
Weight of product 0.0022 kg
Weight of case 9.05 kg
EAN of inner box 4038917126489
EAN case 4038917080347

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

