



Product description

Order number	70.1186.210
Product description	Filter tip, working volume: 1250 µl, transparent, filling level rings, Biosphere®, suitable for 1000 µl pipettes and 1200 µl multi-channel pipettes, 96 piece(s)/box

Product characteristics

Graduation	yes
Suitable for	1000 µl pipettes and 1200 µl multi-channel pipettes
Volume of work	1250 µl
Filter	yes
Coding colour	blue

Size

Max. volume	1250 µl
Length of product	103 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
**Filter tip, 1250 µl, transparent, Biosphere® , 96
piece(s)/box**

Page 2

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Filter material	Polyethylene (PE)

Purity & certification

Product category	no medical device IVD
Purity standard	Biosphere®
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order qty.	768
Type of inner packaging	box
Piece(s) / inner box	96
Piece(s) / case	768
Piece(s) / pallet	90000
Case volume	0,0152 cbm
Weight of product	0,0022 kg
Weight of case	2,08 kg
EAN of inner box	4038917126496
EAN case	4038917080354

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-01-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT