



Product description

Order number	70.1186.410
Product description	Filter tip, working volume: 1250 µl, transparent, filling level rings, PCR Performance Tested, suitable for 1000 µl pipettes and 1200 µl multi-channel pipettes, 384 piece(s)/bag

Product characteristics

Graduation	yes
Suitable for	1000 µl pipettes and 1200 µl multi-channel pipettes
Volume of work	1250 µl
Filter	yes
Coding colour	blue

Size

Max. volume	1250 µl
Length of product	103 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Filter material	Polyethylene (PE)

Purity & certification

Product category	no medical device IVD
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Minimum order qty.	3840
Type of inner packaging	bag
Piece(s) / inner box	384
Piece(s) / case	3840
Piece(s) / pallet	134400
Case volume	0,0374 cbm
Weight of product	0,0011 kg
Weight of case	4,82 kg
EAN of inner box	4038917126502
EAN case	4038917080361

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

