



Product description

Order number	70.1187.101
Product description	Pipette tip, working volume: 10 ml, transparent, filling level rings, suitable for Eppendorf, Gilson and Brand and products of identical design, 25 piece(s)/box

Product characteristics

Graduation	yes
Suitable for	Eppendorf, Gilson and Brand and products of identical design
Volume of work	10 ml
Filter	no
Autoclavable (max. °C)	121 °C

Size

Max. volume	10 ml
Length of product	160.8 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Product category not regulated
CE certified no CE mark
Batched yes

Packaging

Minimum order qty. 250
Type of smallest
subpackaging box
Piece(s) / inner box 25
Piece (s) / outer case 250
Piece(s) / pallet 7500
Depth of box 150 mm
Width of box 118 mm
Height of box 167 mm
Depth of case 610 mm
Width of case 305 mm
Height of case 178 mm
Case volume 0.033 cbm
Weight of product 0.0125 kg
Weight of case 3.63 kg
EAN of inner box 4038917126526
EAN case 4038917080422

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

