



Product description

Order number	70.1187.101
Product description	Pipette tip, working volume: 10 ml, transparent, filling level rings, suitable for Eppendorf, Gilson and Brand and products of identical design, 25 piece(s)/box

Product characteristics

Graduation	yes
Suitable for	Eppendorf, Gilson and Brand and products of identical design
Volume of work	10 ml
Filter	no
Autoclavable (max. °C)	121 °C

Size

Max. volume	10 ml
Length of product	144 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 250
Type of inner packaging box
Piece(s) / inner box 25
Piece(s) / case 250
Piece(s) / pallet 7500
Case volume 0,0316 cbm
Weight of product 0,0125 kg
Weight of case 3,63 kg
EAN of inner box 4038917126526
EAN case 4038917080422

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

