



### **Product description**

Order number	70.1187.101
Product description	Pipette tip, working volume: 10 ml, transparent, filling level rings, suitable for Eppendorf, Gilson and Brand and products of identical design, 25 piece(s)/box

### **Product characteristics**

Graduation	yes
Suitable for	Eppendorf, Gilson and Brand and products of identical design
Volume of work	10 ml
Filter	no
Autoclavable (max. °C)	121 °C

### **Size**

Max. volume	10 ml
Length of product	160.8 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

**Purity & certification**

Product category	no medical device   IVD
Batched	yes

**Packaging**

Minimum order qty.	250
Type of smallest subpackaging	box
Piece(s) / inner box	25
Piece (s) / outer case	250
Piece(s) / pallet	7500
Depth of case	305 mm
Width of case	178 mm
Height of case	610 mm
Case volume	0.033 cbm
Weight of product	0.0125 kg
Weight of case	3.63 kg
EAN of inner box	4038917126526
EAN case	4038917080422

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

