



Product description

Order number	70.1187.102
Product description	Pipette tip, working volume: 10 ml, blue, filling level rings, suitable for Sarpette® M, 25 piece(s)/box

Product characteristics

Graduation	yes
Suitable for	Sarpette® M
Volume of work	10 ml
Filter	no
Autoclavable (max. °C)	121 °C

Size

Max. volume	10 ml
Length of product	144 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Pipette tip, 10 ml, blue, 25 piece(s)/box

Page 2

Material & colours

Product material	Polypropylene (PP)
Colour of product	blue

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	250
Type of smallest subpackaging	box
Piece(s) / inner box	25
Piece (s) / outer case	250
Piece(s) / pallet	7500
Depth of case	603 mm
Width of case	303 mm
Height of case	152 mm
Case volume	0.0278 cbm
Weight of product	0.0128 kg
Weight of case	3.6 kg
EAN of inner box	4038917951494
EAN case	4038917951500

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-27

This document was prepared by EDP support and is valid without signature.



SARSTEDT