




**Product Specification**  
**Gelloader pipette tip, 200 µl, transparent, 96  
piece(s)/box**

**Page 1**



**Product description**

Order number	70.1190.100
Product description	Gelloader pipette tip, working volume: 200 µl, transparent, suitable for common pipettors, 96 piece(s)/box

**Product characteristics**

Graduation	no
Example of use	Gelloader
Suitable for	common pipettors
Volume of work	200 µl
Filter	no
Coding colour	yellow
Autoclavable (max. °C)	121 °C

**Size**

Max. volume	200 µl
Length of product	67 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Gelloader pipette tip, 200 µl, transparent, 96  
piece(s)/box**

**Page 2**

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

**Purity & certification**

Product category	no medical device   IVD
Batched	no

**Packaging**

Minimum order qty.	960
Type of smallest subpackaging	box
Piece(s) / inner box	96
Piece (s) / outer case	960
Piece(s) / pallet	103680
Depth of case	475 mm
Width of case	255 mm
Height of case	100 mm
Case volume	0.0121 cbm
Weight of product	0.0015 kg
Weight of case	1.73 kg
EAN of inner box	4038917344593
EAN case	4038917344586

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-08

This document was prepared by EDP  
support and is valid without signature.



**SARSTEDT**