



### **Product description**

Order number	70.1193.150
Product description	Filter tip, working volume: 1000 µl, black, suitable for application on SARSTEDT IVARO AP devices, 480 piece(s)/blister

### **Product characteristics**

Suitable for	application on SARSTEDT IVARO AP devices
Volume of work	1000 µl
Filter	yes

### **Size**

Max. volume	1000 µl
Length of product	95 mm

### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	black

### **Purity & certification**

Product category	IVD product
CE certified	CE supplier > no Sarstedt CE
Batched	yes

### **Packaging**

Minimum order qty.	3840
Type of inner packaging	blister
Piece(s) / inner box	3840
Piece(s) / case	3840
Piece(s) / pallet	61440
Case volume	0,0754 cbm
Weight of product	0,0013 kg
Weight of case	7,9 kg
EAN of inner box	4038917434768
EAN case	4038917434751

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

