



Product description

| | |
|---------------------|--|
| Order number | 70.1193.350 |
| Product description | Filter tip, working volume: 300 µl, black, suitable for application on SARSTEDT IVARO AP devices, 480 piece(s)/blister |

Product characteristics

| | |
|----------------|--|
| Suitable for | application on SARSTEDT IVARO AP devices |
| Volume of work | 300 µl |
| Filter | yes |

Size

| | |
|-------------|--------|
| Max. volume | 300 µl |
|-------------|--------|

Material & colours

| | |
|-------------------|--------------------|
| Product material | Polypropylene (PP) |
| Colour of product | black |

Purity & certification

| | |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified | CE supplier |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 3840 |
| Type of smallest subpackaging | blister |
| Piece(s) / inner box | 480 |
| Piece (s) / outer case | 3840 |
| Piece(s) / pallet | 61440 |
| Depth of case | 488 mm |
| Width of case | 430 mm |
| Height of case | 325 mm |
| Case volume | 0.0682 cbm |
| Weight of product | 0.0018 kg |
| Weight of case | 6.87 kg |
| EAN of inner box | 4038917898539 |
| EAN case | 4038917898546 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

