



Product description

Order number	70.1193.350
Product description	Filter tip, working volume: 300 µl, black, suitable for application on SARSTEDT IVARO AP devices, 480 piece(s)/blister

Product characteristics

Suitable for	application on SARSTEDT IVARO AP devices
Volume of work	300 µl
Filter	yes

Size

Max. volume	300 µl
-------------	--------

Material & colours

Product material	Polypropylene (PP)
Colour of product	black

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE supplier
Batched	yes

Packaging

Minimum order qty.	3840
Type of smallest subpackaging	blister
Piece(s) / inner box	480
Piece (s) / outer case	3840
Piece(s) / pallet	61440
Depth of case	488 mm
Width of case	430 mm
Height of case	325 mm
Case volume	0.0682 cbm
Weight of product	0.0018 kg
Weight of case	6.87 kg
EAN of inner box	4038917898539
EAN case	4038917898546

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

