



Product description

| | |
|---------------------|--|
| Order number | 70.3010.285 |
| Product description | Filter tip, working volume: 2.5 µl, transparent, filling level rings, Biosphere® plus, Low retention, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpiquette and Brand and products of identical design, 96 piece(s)/box |

Product characteristics

| | |
|----------------|--|
| Graduation | yes |
| Suitable for | SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpiquette and Brand and products of identical design |
| Volume of work | 2,5 µl |
| Filter | yes |
| Coding colour | grey |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Size

| | |
|-------------------|--------|
| Max. volume | 2.5 µl |
| Length of product | 31 mm |

Material & colours

| | |
|-------------------|--------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |
| Filter material | Polyethylene (PE) |

Purity & certification

| | |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified | CE |
| Purity standard | Biosphere® plus |
| Sterilisation | Gas sterilisation |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1920 |
| Type of smallest subpackaging | box |
| Piece(s) / inner box | 480 |
| Piece (s) / outer case | 1920 |
| Piece(s) / pallet | 96000 |
| Depth of box | 330 mm |
| Width of box | 111 mm |
| Height of box | 148 mm |
| Depth of case | 348 mm |
| Width of case | 325 mm |
| Height of case | 245 mm |
| Case volume | 0.0277 cbm |
| Weight of product | 0.0009 kg |
| Weight of case | 5.16 kg |
| EAN of inner box | 4038917894661 |
| EAN case | 4038917894678 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

