



Product description

| | |
|---------------------|--|
| Order number | 70.3010.365 |
| Product description | Filter tip, working volume: 2.5 µl, transparent, filling level rings, Biosphere® plus, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpiette and Brand and products of identical design, 96 piece(s)/SingleRefill |

Product characteristics

| | |
|----------------|--|
| Graduation | yes |
| Suitable for | SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpiette and Brand and products of identical design |
| Volume of work | 2,5 µl |
| Filter | yes |
| Coding colour | grey |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Size

| | |
|-------------------|--------|
| Max. volume | 2.5 µl |
| Length of product | 31 mm |

Material & colours

| | |
|-------------------|--------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |
| Filter material | Polyethylene (PE) |

Purity & certification

| | |
|------------------|--------------------------------|
| Product category | In-vitro Diagnostikum (IVD) |
| CE certified | CE - HerstellerSelbsterklärung |
| Purity standard | Biosphere® plus |
| Sterilisation | Gas sterilisation |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1920 |
| Type of smallest subpackaging | SingleRefill |
| Piece(s) / inner box | 480 |
| Piece (s) / outer case | 1920 |
| Piece(s) / pallet | 115200 |
| Depth of box | 530 mm |
| Width of box | 111 mm |
| Height of box | 148 mm |
| Depth of case | 437 mm |
| Width of case | 307 mm |
| Height of case | 170 mm |
| Case volume | 0.0436 cbm |
| Weight of product | 0.0006 kg |
| Weight of case | 2.18 kg |
| EAN of inner box | 4038917895040 |
| EAN case | 4038917895057 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

