



### Product description

Order number	70.3010
Product description	Pipettenspitze, Arbeitsvolumen: 10 µl, transparent, Füllstandsringe, PCR Performance Tested, passend für Eppendorf, Gilson, Finnpiquette, Brand und Socorex sowie baugleiche Ausführungen, 1000 Stück/Beutel

### Product characteristics

Graduation	yes
Suitable for	Eppendorf, Gilson, Finnpiquette, Brand and Socorex and products of identical design
Volume of work	10 µl
Filter	no
Coding colour	grey
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Product Specification**  
**Pipette tip, 10 µl, transparent, PCR Performance**  
**Tested, 1000 piece(s)/bag**

**Page 2**

**Size**

Max. volume 10 µl  
Length of product 31 mm

**Material & colours**

Product material Polypropylene (PP)  
Colour of product transparent

**Purity & certification**

Product in accordance with IVD product  
MDD or IVDD  
CE certified CE  
Purity standard PCR Performance Tested  
Lot or serial number yes

**Packaging**

Minimum order qty. 10000  
Type of inner packaging bag  
Piece(s) / inner box 1000  
Piece(s) / case 10000  
Piece(s) / pallet 960000  
Depth of case 348 mm  
Width of case 278 mm  
Height of case 140 mm  
Case volume 0,0135 cbm  
Weight of product 0,0001 kg  
Weight of case 1,64 kg  
EAN of inner box 4038917525824  
EAN case 4038917525817

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2021-05-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**