



**Product description**

|                     |                                                                                                                                                                                                                                    |
|---------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Order number        | 70.3020.100                                                                                                                                                                                                                        |
| Product description | Pipette tip, working volume: 20 µl, transparent, filling level rings, PCR Performance Tested, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette and Brand and products of identical design, 480 piece(s)/StackPack |

**Product characteristics**

|                |                                                                                                 |
|----------------|-------------------------------------------------------------------------------------------------|
| Graduation     | yes                                                                                             |
| Suitable for   | SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette and Brand and products of identical design |
| Volume of work | 20 µl                                                                                           |
| Filter         | no                                                                                              |
| Coding colour  | grey                                                                                            |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Pipette tip, 20 µl, transparent, PCR Performance**  
**Tested, 480 piece(s)/StackPack**

**Page 2**

**Size**

|                   |       |
|-------------------|-------|
| Max. volume       | 20 µl |
| Length of product | 46 mm |

**Material & colours**

|                   |                    |
|-------------------|--------------------|
| Product material  | Polypropylene (PP) |
| Colour of product | transparent        |

**Purity & certification**

|                  |                               |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified     | CE                            |
| Purity standard  | PCR Performance Tested        |
| Batched          | yes                           |

**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 1920          |
| Type of smallest subpackaging | StackPack     |
| Piece(s) / inner box          | 1920          |
| Piece (s) / outer case        | 1920          |
| Piece(s) / pallet             | 215040        |
| Depth of case                 | 290 mm        |
| Width of case                 | 216 mm        |
| Height of case                | 225 mm        |
| Case volume                   | 0.0141 cbm    |
| Weight of product             | 0.0005 kg     |
| Weight of case                | 1.69 kg       |
| EAN of inner box              | 4038917893978 |
| EAN case                      | 4038917893985 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-06

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**