



### Product description

Order number	70.3020
Product description	Pipette tip, working volume: 20 µl, transparent, filling level rings, PCR Performance Tested, suitable for Eppendorf, Gilson, Finnpiquette, Brand and Socorex and products of identical design, 1000 piece(s)/bag

### Product characteristics

Graduation	yes
Suitable for	Eppendorf, Gilson, Finnpiquette, Brand and Socorex and products of identical design
Volume of work	20 µl
Filter	no
Coding colour	grey
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### Size

Max. volume	20 µl
Length of product	46 mm

### Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

### Purity & certification

Product category	IVD product
CE certified	CE
Purity standard	PCR Performance Tested
Batched	yes

### Packaging

Minimum order qty.	10000
Type of inner packaging	bag
Piece(s) / inner box	1000
Piece(s) / case	10000
Piece(s) / pallet	720000
Depth of case	533 mm
Width of case	263 mm
Height of case	140 mm
Case volume	0,0196 cbm
Weight of product	0,0002 kg
Weight of case	2,19 kg
EAN of inner box	4038917525862
EAN case	4038917525855

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

