



### Product description

Order number	70.3030.200
Product description	Pipette tip, working volume: 200 µl, transparent, filling level rings, PCR Performance Tested, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design, 96 piece(s)/box

### Product characteristics

Graduation	yes
Suitable for	SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design
Volume of work	200 µl
Filter	no
Coding colour	yellow
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Size**

Max. volume	200 µl
Length of product	51 mm

### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

### **Purity & certification**

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	PCR Performance Tested
Batched	yes

### **Packaging**

Minimum order qty.	1920
Type of smallest subpackaging	box
Piece(s) / inner box	480
Piece (s) / outer case	1920
Piece(s) / pallet	96000
Depth of box	330 mm
Width of box	111 mm
Height of box	148 mm
Depth of case	348 mm
Width of case	325 mm
Height of case	245 mm
Case volume	0.0277 cbm
Weight of product	0.0011 kg
Weight of case	3.12 kg
EAN of inner box	4038917847117
EAN case	4038917847124

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

