



### Product description

|                     |                                                                                                                                                                                                                                      |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Order number        | 70.3030.200                                                                                                                                                                                                                          |
| Product description | Pipette tip, working volume: 200 µl, transparent, filling level rings, PCR Performance Tested, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design, 96 piece(s)/box |

### Product characteristics

|                        |                                                                                                         |
|------------------------|---------------------------------------------------------------------------------------------------------|
| Graduation             | yes                                                                                                     |
| Suitable for           | SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design |
| Volume of work         | 200 µl                                                                                                  |
| Filter                 | no                                                                                                      |
| Coding colour          | yellow                                                                                                  |
| Autoclavable (max. °C) | 121 °C                                                                                                  |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Size**

|                   |        |
|-------------------|--------|
| Max. volume       | 200 µl |
| Length of product | 51 mm  |

### **Material & colours**

|                   |                    |
|-------------------|--------------------|
| Product material  | Polypropylene (PP) |
| Colour of product | transparent        |

### **Purity & certification**

|                  |                                      |
|------------------|--------------------------------------|
| Product category | In-vitro diagnostic (IVD)            |
| CE certified     | CE - manufacturer's self-declaration |
| Purity standard  | PCR Performance Tested               |
| Batched          | yes                                  |

### **Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 1920          |
| Type of smallest subpackaging | box           |
| Piece(s) / inner box          | 480           |
| Piece (s) / outer case        | 1920          |
| Piece(s) / pallet             | 96000         |
| Depth of box                  | 330 mm        |
| Width of box                  | 111 mm        |
| Height of box                 | 148 mm        |
| Depth of case                 | 348 mm        |
| Width of case                 | 325 mm        |
| Height of case                | 245 mm        |
| Case volume                   | 0.0277 cbm    |
| Weight of product             | 0.0011 kg     |
| Weight of case                | 3.12 kg       |
| EAN of inner box              | 4038917847117 |
| EAN case                      | 4038917847124 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

