



Product description

Order number	70.3030.305
Product description	Pipette tip, working volume: 200 µl, transparent, filling level rings, Biosphere® plus, suitable for Eppendorf, Gilson, Finnpiquette, Biohit, Brand and Socorex and products of identical design, 96 piece(s)/SingleRefill

Product characteristics

Graduation	yes
Suitable for	Eppendorf, Gilson, Finnpiquette, Biohit, Brand and Socorex and products of identical design
Volume of work	200 µl
Filter	no
Coding colour	yellow
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Size

Max. volume	200 µl
Length of product	51 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product category	In-vitro diagnostic, CE
CE certified	CE
Purity standard	Biosphere® plus
Sterilisation	Gas sterilisation
Batched	yes

Packaging

Minimum order qty.	1920
Type of inner packaging	SingleRefill
Piece(s) / inner box	480
Piece(s) / case	1920
Piece(s) / pallet	115200
Depth of case	437 mm
Width of case	307 mm
Height of case	170 mm
Case volume	0,0228 cbm
Weight of product	0,0008 kg
Weight of case	5,16 kg
EAN of inner box	4038917847032
EAN case	4038917847049

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

