




Product Specification
Filter tip, 20 µl, transparent, Biosphere® plus,
Low retention, 96 piece(s)/SingleRefill

Page 1



Product description

| | |
|---------------------|---|
| Order number | 70.3030.385 |
| Product description | Filter tip, working volume: 20 µl, transparent, filling level rings, Biosphere® plus, Low retention, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design, 96 piece(s)/SingleRefill |

Product characteristics

| | |
|----------------|---|
| Graduation | yes |
| Suitable for | SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design |
| Volume of work | 20 µl |
| Filter | yes |
| Coding colour | yellow |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

**Filter tip, 20 µl, transparent, Biosphere® plus,
Low retention, 96 piece(s)/SingleRefill**

Page 2

Size

| | |
|-------------------|-------|
| Max. volume | 20 µl |
| Length of product | 51 mm |

Material & colours

| | |
|-------------------|--------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |
| Filter material | Polyethylene (PE) |

Purity & certification

| | |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified | CE |
| Purity standard | Biosphere® plus |
| Sterilisation | Gas sterilisation |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1920 |
| Type of smallest subpackaging | SingleRefill |
| Piece(s) / inner box | 480 |
| Piece (s) / outer case | 1920 |
| Piece(s) / pallet | 115200 |
| Depth of box | 291 mm |
| Width of box | 105 mm |
| Height of box | 148 mm |
| Depth of case | 437 mm |
| Width of case | 307 mm |
| Height of case | 170 mm |
| Case volume | 0.0228 cbm |
| Weight of product | 0.0008 kg |
| Weight of case | 5.16 kg |
| EAN of inner box | 4038917847070 |
| EAN case | 4038917847087 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-03

This document was prepared by EDP support and is valid without signature.



SARSTEDT