



### Product description

Order number	70.3030
Product description	Pipette tip, working volume: 200 µl, transparent, filling level rings, PCR Performance Tested, suitable for Eppendorf, Gilson, Finnpipette, Biohit, Brand and Socorex and products of identical design, 1000 piece(s)/bag

### Product characteristics

Graduation	yes
Suitable for	Eppendorf, Gilson, Finnpipette, Biohit, Brand and Socorex and products of identical design
Volume of work	200 µl
Filter	no
Coding colour	yellow
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### Size

Max. volume	200 µl
Length of product	51 mm

### Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

### Purity & certification

Product category	IVD product
CE certified	CE
Purity standard	PCR Performance Tested
Batched	yes

### Packaging

Minimum order qty.	10000
Type of inner packaging	bag
Piece(s) / inner box	1000
Piece(s) / case	10000
Piece(s) / pallet	480000
Depth of case	388 mm
Width of case	378 mm
Height of case	215 mm
Case volume	0,0315 cbm
Weight of product	0,0003 kg
Weight of case	3,83 kg
EAN of inner box	4038917525947
EAN case	4038917525930

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

