



## Product description

Order number	70.3031.210
Product description	Pipette tip, working volume: 250 µl, transparent, filling level rings, PCR Performance Tested, Low retention, suitable for Eppendorf, Gilson, Finnpipette, Biohit, Brand and Socorex and products of identical design, 96 piece(s)/box

## Product characteristics

Graduation	yes
Suitable for	Eppendorf, Gilson, Finnpipette, Biohit, Brand and Socorex and products of identical design
Volume of work	250 µl
Filter	no
Coding colour	yellow
Autoclavable (max. °C)	121 °C

## Size

Max. volume	250 µl
Length of product	51 mm

## Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

## Purity & certification

Product category	IVD product
CE certified	CE
Purity standard	PCR Performance Tested
Batched	yes

## Packaging

Minimum order qty.	1920
Type of inner packaging	box
Piece(s) / inner box	480
Piece(s) / case	1920
Piece(s) / pallet	96000
Depth of case	348 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Product Specification**  
**Pipette tip, 250 µl, transparent, PCR Performance**  
**Tested, Low retention, 96 piece(s)/box**

**Page 2**

Width of case	325 mm
Height of case	245 mm
Case volume	0,0277 cbm
Weight of product	0,0011 kg
Weight of case	5,16 kg
EAN of inner box	4038917845762
EAN case	4038917845779

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2021-12-03

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**