



### Product description

Order number	70.3031.305
Product description	Pipette tip, working volume: 250 µl, transparent, filling level rings, Biosphere® plus, suitable for Eppendorf, Gilson, Finnpiquette, Biohit, Brand and Socorex and products of identical design, 96 piece(s)/SingleRefill

### Product characteristics

Graduation	yes
Suitable for	Eppendorf, Gilson, Finnpiquette, Biohit, Brand and Socorex and products of identical design
Volume of work	250 µl
Filter	no
Coding colour	yellow
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Product Specification**  
**Pipette tip, 250 µl, transparent, Biosphere® plus,**  
**96 piece(s)/SingleRefill**

**Page 2**

**Size**

Max. volume 250 µl  
Length of product 51 mm

**Material & colours**

Product material Polypropylene (PP)  
Colour of product transparent

**Purity & certification**

Product category IVD product  
CE certified CE  
Purity standard Biosphere® plus  
Sterilisation Gas sterilisation  
Batched yes

**Packaging**

Minimum order qty. 1920  
Type of inner packaging SingleRefill  
Piece(s) / inner box 480  
Piece(s) / case 1920  
Piece(s) / pallet 115200  
Depth of case 437 mm  
Width of case 307 mm  
Height of case 170 mm  
Case volume 0,0228 cbm  
Weight of product 0,0008 kg  
Weight of case 5,16 kg  
EAN of inner box 4038917845922  
EAN case 4038917845939

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2021-12-03

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**