



Product description

Order number	70.3040.100
Product description	Pipette tip, working volume: 300 µl, transparent, filling level rings, PCR Performance Tested, suitable for Eppendorf, Gilson, Finnpipette, Biohit, Brand and Socorex and products of identical design, 480 piece(s)/StackPack

Product characteristics

Graduation	yes
Suitable for	Eppendorf, Gilson, Finnpipette, Biohit, Brand and Socorex and products of identical design
Volume of work	300 µl
Filter	no
Coding colour	orange
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Size

Max. volume	300 µl
Diameter	7 mm
Length of product	60 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product category	IVD product
Certificate of conformity	yes
CE certified	CE
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Minimum order qty.	1920
Type of inner packaging	StackPack
Piece(s) / inner box	1920
Piece(s) / case	1920
Piece(s) / pallet	245760
Depth of case	290 mm
Width of case	186 mm
Height of case	225 mm
Case volume	0,0121 cbm
Weight of product	0,0007 kg
Weight of case	1,69 kg
EAN of inner box	4038917433433
EAN case	4038917433426

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

