



**Product description**

Order number	70.3040.110
Product description	Pipette tip, working volume: 300 µl, transparent, filling level rings, PCR Performance Tested, Low retention, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design, 480 piece(s)/StackPack

**Product characteristics**

Graduation	yes
Suitable for	SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design
Volume of work	300 µl
Filter	no
Coding colour	orange

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

**Pipette tip, 300 µl, transparent, PCR Performance Tested, Low retention, 480 piece(s)/StackPack**

**Page 2**

### Size

Max. volume	300 µl
Diameter	7 mm
Length of product	60 mm

### Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

### Purity & certification

Product category	In vitrodiagnostic device, CE
Certificate of conformity	yes
CE certified	CE
Purity standard	PCR Performance Tested
Batched	yes

### Packaging

Minimum order qty.	1920
Type of smallest subpackaging	StackPack
Piece(s) / inner box	1920
Piece (s) / outer case	1920
Piece(s) / pallet	245760
Depth of case	290 mm
Width of case	186 mm
Height of case	225 mm
Case volume	0.0121 cbm
Weight of product	0.0007 kg
Weight of case	1.69 kg
EAN of inner box	4038917433402
EAN case	4038917433396

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-05

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**