



### **Product description**

Order number	70.3040.200
Product description	Pipette tip, working volume: 300 µl, transparent, filling level rings, PCR Performance Tested, suitable for Eppendorf, Gilson, Finnpiquette, Biohit and Brand and products of identical design, 96 piece(s)/box

### **Product characteristics**

Graduation	yes
Suitable for	Eppendorf, Gilson, Finnpiquette, Biohit and Brand and products of identical design
Volume of work	300 µl
Filter	no
Coding colour	orange
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### Size

Max. volume	300 µl
Diameter	7 mm
Length of product	60 mm

### Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

### Purity & certification

Product category	In-vitro diagnostic, CE
Certificate of conformity	yes
CE certified	CE
Purity standard	PCR Performance Tested
Batched	yes

### Packaging

Minimum order qty.	1920
Type of inner packaging	box
Piece(s) / inner box	480
Piece(s) / case	1920
Piece(s) / pallet	61440
Depth of case	547 mm
Width of case	325 mm
Height of case	245 mm
Case volume	0,0436 cbm
Weight of product	0,0013 kg
Weight of case	4,05 kg
EAN of inner box	4038917433280
EAN case	4038917433273

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

