



### Product description

Order number	70.3040.355
Product description	Filter tip, working volume: 300 µl, transparent, filling level rings, Biosphere® plus, suitable for Eppendorf, Gilson and Finnpiquette and products of identical design, 96 piece(s)/SingleRefill

### Product characteristics

Graduation	yes
Suitable for	Eppendorf, Gilson and Finnpiquette and products of identical design
Volume of work	300 µl
Filter	yes
Coding colour	orange

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### Size

Max. volume	300 µl
Diameter	7 mm
Length of product	60 mm

### Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

### Purity & certification

Product category	IVD product
Certificate of conformity	yes
CE certified	CE
Purity standard	Biosphere® plus
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Batched	yes

### Packaging

Minimum order qty.	1920
Type of inner packaging	SingleRefill
Piece(s) / inner box	480
Piece(s) / case	1920
Piece(s) / pallet	103680
Depth of case	437 mm
Width of case	351 mm
Height of case	172 mm
Case volume	0,0264 cbm
Weight of product	0,0009 kg
Weight of case	2,81 kg
EAN of inner box	4038917433198
EAN case	4038917433181

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

