



### Product description

|                     |   |
|---------------------|---|
| Order number        | 70.3040.375   |
| Product description | Filter tip, working volume: 300 µl, transparent, filling level rings, Biosphere® plus, Low retention, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson and Brand and products of identical design, 96 piece(s)/SingleRefill |

### Product characteristics

|                |  |
|----------------|--|
| Graduation     | yes  |
| Suitable for   | SARSTEDT Sarpette® M, Eppendorf, Gilson and Brand and products of identical design |
| Volume of work | 300 µl   |
| Filter         | yes  |
| Coding colour  | orange   |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

**Filter tip, 300 µl, transparent, Biosphere® plus,  
Low retention, 96 piece(s)/SingleRefill**

**Page 2**

### Size

|                   |        |
|-------------------|--------|
| Max. volume       | 300 µl |
| Diameter          | 7 mm   |
| Length of product | 60 mm  |

### Material & colours

|                   |                    |
|-------------------|--------------------|
| Product material  | Polypropylene (PP) |
| Colour of product | transparent        |

### Purity & certification

|                           |                                |
|---------------------------|--------------------------------|
| Product category          | In vitro diagnostic device, CE |
| Certificate of conformity | yes                            |
| CE certified              | CE                             |
| Purity standard           | Biosphere® plus                |
| Sterilisation             | Gas sterilisation              |
| Pyrogen-free              | yes                            |
| Batched                   | yes                            |

### Packaging

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 1920          |
| Type of smallest subpackaging | SingleRefill  |
| Piece(s) / inner box          | 480           |
| Piece (s) / outer case        | 1920          |
| Piece(s) / pallet             | 103680        |
| Depth of box                  | 333 mm        |
| Width of box                  | 105 mm        |
| Height of box                 | 148 mm        |
| Depth of case                 | 437 mm        |
| Width of case                 | 351 mm        |
| Height of case                | 172 mm        |
| Case volume                   | 0.0264 cbm    |
| Weight of product             | 0.0009 kg     |
| Weight of case                | 2.81 kg       |
| EAN of inner box              | 4038917433167 |
| EAN case                      | 4038917433150 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-29

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**