



### **Product description**

Order number	70.3050.020
Product description	Pipette tip, working volume: 1000 µl, blue, filling level rings, PCR Performance Tested, suitable for Eppendorf, Gilson, Finnpiquette, Biohit, Brand and Socorex and products of identical design, 500 piece(s)/bag

### **Product characteristics**

Graduation	yes
Suitable for	Eppendorf, Gilson, Finnpiquette, Biohit, Brand and Socorex and products of identical design
Volume of work	1000 µl
Filter	no
Coding colour	blue
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Product Specification**  
**Pipette tip, 1000 µl, blue, PCR Performance**  
**Tested, 500 piece(s)/bag**

**Page 2**

**Size**

Max. volume	1000 µl
Diameter	9 mm
Length of product	72 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	blue

**Purity & certification**

Product category	IVD product
Certificate of conformity	yes
CE certified	CE
Purity standard	PCR Performance Tested
Batched	yes

**Packaging**

Minimum order qty.	5000
Type of inner packaging	bag
Piece(s) / inner box	500
Piece(s) / case	5000
Piece(s) / pallet	180000
Depth of case	398 mm
Width of case	398 mm
Height of case	300 mm
Case volume	0,0475 cbm
Weight of product	0,0008 kg
Weight of case	4,77 kg
EAN of inner box	4038917434584
EAN case	4038917434577

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2021-12-03

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**