




Product Specification
Pipette tip, 1000 µl, black, PCR Performance
Tested, 500 piece(s)/bag

Page 1



Product description

Order number	70.3050.030
Product description	Pipette tip, working volume: 1000 µl, black, PCR Performance Tested, 500 piece(s)/bag

Product characteristics

Graduation	no
Volume of work	1000 µl
Filter	no
Coding colour	blue
Autoclavable (max. °C)	121 °C

Size

Max. volume	1000 µl
Diameter	9 mm
Length of product	72 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-01-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
Pipette tip, 1000 µl, black, PCR Performance
Tested, 500 piece(s)/bag

Page 2

Material & colours

Product material Polypropylene (PP)
Colour of product black

Purity & certification

Product category no medical device | IVD
Certificate of conformity yes
Purity standard PCR Performance Tested
Batched yes

Packaging

Minimum order qty. 5000
Type of inner packaging bag
Piece(s) / inner box 500
Piece(s) / case 5000
Piece(s) / pallet 180000
Depth of case 398 mm
Width of case 398 mm
Height of case 300 mm
Case volume 0,0475 cbm
Weight of product 0,0009 kg
Weight of case 5,24 kg
EAN of inner box 4038917434614
EAN case 4038917434607

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-01-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT