




**Product Specification**  
**Pipette tip, 1000 µl, transparent, PCR**  
**Performance Tested, Low retention, 480**  
**piece(s)/StackPack**

**Page 1**



**Product description**

Order number	70.3050.110
Product description	Pipette tip, working volume: 1000 µl, transparent, filling level rings, PCR Performance Tested, Low retention, suitable for Eppendorf, Gilson, Finnpiquette, Biohit, Brand and Socorex and products of identical design, 480 piece(s)/StackPack

**Product characteristics**

Graduation	yes
Suitable for	Eppendorf, Gilson, Finnpiquette, Biohit, Brand and Socorex and products of identical design
Volume of work	1000 µl
Filter	no
Coding colour	blue
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2021-12-04

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**Pipette tip, 1000 µl, transparent, PCR**  
**Performance Tested, Low retention, 480**  
**piece(s)/StackPack**

**Page 2**

**Size**

Max. volume	1000 µl
Diameter	9 mm
Length of product	72 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

**Purity & certification**

Product category	IVD product
Certificate of conformity	yes
CE certified	CE
Purity standard	PCR Performance Tested
Batched	yes

**Packaging**

Minimum order qty.	1920
Type of inner packaging	StackPack
Piece(s) / inner box	1920
Piece(s) / case	1920
Piece(s) / pallet	149760
Depth of case	303 mm
Width of case	295 mm
Height of case	225 mm
Case volume	0,0201 cbm
Weight of product	0,0012 kg
Weight of case	2,79 kg
EAN of inner box	4038917424813
EAN case	4038917424806

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2021-12-04

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**