



### Product description

Order number	70.3050.205
Product description	Pipette tip, working volume: 1000 µl, transparent, filling level rings, Biosphere® plus, suitable for Eppendorf, Gilson, Finnpiquette and Brand and products of identical design, 96 piece(s)/box

### Product characteristics

Graduation	yes
Suitable for	Eppendorf, Gilson, Finnpiquette and Brand and products of identical design
Volume of work	1000 µl
Filter	no
Coding colour	blue
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Product Specification**  
**Pipette tip, 1000 µl, transparent, Biosphere® plus,**  
**96 piece(s)/box**

**Page 2**

**Size**

Max. volume	1000 µl
Diameter	9 mm
Length of product	72 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

**Purity & certification**

Product category	IVD product
Certificate of conformity	yes
CE certified	CE
Purity standard	Biosphere® plus
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Batched	yes

**Packaging**

Minimum order qty.	1920
Type of inner packaging	box
Piece(s) / inner box	480
Piece(s) / case	1920
Piece(s) / pallet	61440
Depth of case	547 mm
Width of case	325 mm
Height of case	245 mm
Case volume	0,0436 cbm
Weight of product	0,0018 kg
Weight of case	4,99 kg
EAN of inner box	4038917424752
EAN case	4038917424745

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2021-12-03

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**