



Product description

Order number	70.3050.305
Product description	Pipette tip, working volume: 1000 µl, transparent, filling level rings, Biosphere® plus, suitable for Eppendorf, Gilson, Finnpipette and Brand and products of identical design, 96 piece(s)/SingleRefill

Product characteristics

Graduation	yes
Suitable for	Eppendorf, Gilson, Finnpipette and Brand and products of identical design
Volume of work	1000 µl
Filter	no
Coding colour	blue
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
**Pipette tip, 1000 µl, transparent, Biosphere® plus,
96 piece(s)/SingleRefill**

Page 2

Size

Max. volume	1000 µl
Diameter	9 mm
Length of product	72 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product category	In-vitro diagnostic, CE
Certificate of conformity	yes
CE certified	CE
Purity standard	Biosphere® plus
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order qty.	1920
Type of inner packaging	SingleRefill
Piece(s) / inner box	480
Piece(s) / case	1920
Piece(s) / pallet	69120
Depth of case	512 mm
Width of case	325 mm
Height of case	237 mm
Case volume	0,0394 cbm
Weight of product	0,0014 kg
Weight of case	4,12 kg
EAN of inner box	4038917424639
EAN case	4038917424622

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-01-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT