




**Product Specification**  
**Pipette tip, 1,000 µl, transparent, PCR**  
**Performance Tested, 500 piece(s)/bag**

**Page 1**



**Product description**

Order number	70.3050
Product description	Pipette tip, working volume: 1,000 µl, transparent, filling level rings, PCR Performance Tested, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design, 500 piece(s)/bag

**Product characteristics**

Graduation	yes
Suitable for	SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design
Volume of work	1.000 µl
Filter	no
Coding colour	blue
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Pipette tip, 1,000 µl, transparent, PCR**  
**Performance Tested, 500 piece(s)/bag**

**Page 2**

**Size**

Max. volume	1000 µl
Diameter	9 mm
Length of product	72 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

**Purity & certification**

Product category	In vitrodiagnostic device, CE
Certificate of conformity	yes
CE certified	CE
Purity standard	PCR Performance Tested
Batched	yes

**Packaging**

Minimum order qty.	5000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	5000
Piece(s) / pallet	180000
Depth of case	398 mm
Width of case	398 mm
Height of case	300 mm
Case volume	0.0475 cbm
Weight of product	0.0008 kg
Weight of case	4.76 kg
EAN of inner box	4038917434553
EAN case	4038917434546

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-29

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**