



### **Product description**

Order number	70.3060.100
Product description	Pipette tip, working volume: 1250 µl, transparent, filling level rings, PCR Performance Tested, suitable for Eppendorf, Gilson, Finnpipette, Biohit, Brand and Socorex and products of identical design, 480 piece(s)/StackPack

### **Product characteristics**

Graduation	yes
Suitable for	Eppendorf, Gilson, Finnpipette, Biohit, Brand and Socorex and products of identical design
Volume of work	1250 µl
Filter	no
Coding colour	blue
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Product Specification**  
**Pipette tip, 1250 µl, transparent, PCR**  
**Performance Tested, 480 piece(s)/StackPack**

**Page 2**

**Size**

Max. volume	1250 µl
Diameter	9 mm
Length of product	90 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

**Purity & certification**

Product category	IVD product
Certificate of conformity	yes
CE certified	CE
Purity standard	PCR Performance Tested
Batched	yes

**Packaging**

Minimum order qty.	1920
Type of inner packaging	StackPack
Piece(s) / inner box	1920
Piece(s) / case	1920
Piece(s) / pallet	184320
Depth of case	290 mm
Width of case	267 mm
Height of case	225 mm
Case volume	0,0174 cbm
Weight of product	0,0013 kg
Weight of case	2,97 kg
EAN of inner box	4038917432207
EAN case	4038917432191

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2021-12-04

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**