



### Product description

Order number	70.3060.255
Product description	Filter tip, working volume: 1000 µl, transparent, filling level rings, Biosphere® plus, suitable for Eppendorf, Gilson and Brand and products of identical design, 96 piece(s)/box

### Product characteristics

Graduation	yes
Suitable for	Eppendorf, Gilson and Brand and products of identical design
Volume of work	1000 µl
Filter	yes
Coding colour	blue

### Size

Max. volume	1000 µl
Diameter	9 mm
Length of product	90 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Material & colours**

Product material Polypropylene (PP)  
Colour of product transparent

### **Purity & certification**

Product category IVD product  
Certificate of conformity yes  
CE certified CE  
Purity standard Biosphere® plus  
Sterilisation Gas sterilisation  
Pyrogen-free yes  
Batched yes

### **Packaging**

Minimum order qty. 1920  
Type of inner packaging box  
Piece(s) / inner box 480  
Piece(s) / case 1920  
Piece(s) / pallet 61440  
Depth of case 547 mm  
Width of case 325 mm  
Height of case 245 mm  
Case volume 0,0436 cbm  
Weight of product 0,002 kg  
Weight of case 5,41 kg  
EAN of inner box 4038917432269  
EAN case 4038917432252

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

