



Product description

Order number	70.3060.255
Product description	Filter tip, working volume: 1000 µl, transparent, filling level rings, Biosphere® plus, suitable for Eppendorf, Gilson and Brand and products of identical design, 96 piece(s)/box

Product characteristics

Graduation	yes
Suitable for	Eppendorf, Gilson and Brand and products of identical design
Volume of work	1000 µl
Filter	yes
Coding colour	blue

Size

Max. volume	1000 µl
Diameter	9 mm
Length of product	90 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Product category IVD product
Certificate of conformity yes
CE certified CE
Purity standard Biosphere® plus
Sterilisation Gas sterilisation
Pyrogen-free yes
Batched yes

Packaging

Minimum order qty. 1920
Type of inner packaging box
Piece(s) / inner box 480
Piece(s) / case 1920
Piece(s) / pallet 61440
Depth of case 547 mm
Width of case 325 mm
Height of case 245 mm
Case volume 0,0436 cbm
Weight of product 0,002 kg
Weight of case 5,41 kg
EAN of inner box 4038917432269
EAN case 4038917432252

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

