



### **Product description**

Order number	70.3070.205
Product description	Pipette tip, working volume: 1,250 µl, transparent, filling level rings, Biosphere® plus, 96 piece(s)/box

### **Product characteristics**

Graduation	yes
Volume of work	1.250 µl
Filter	no
Coding colour	blue
Autoclavable (max. °C)	121 °C

### **Size**

Max. volume	1250 µl
Diameter	9 mm
Length of product	108 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Pipette tip, 1,250 µl, transparent, Biosphere® plus, 96 piece(s)/box**

**Page 2**

**Material & colours**

Product material Polypropylene (PP)  
Colour of product transparent

**Purity & certification**

Product category In-vitro diagnostic (IVD)  
CE certified CE - manufacturer's self-declaration  
Purity standard Biosphere® plus  
Sterilisation Gas sterilisation  
Batched yes

**Packaging**

Minimum order qty. 1920  
Type of smallest box  
subpackaging  
Piece(s) / inner box 480  
Piece (s) / outer case 1920  
Piece(s) / pallet 46080  
Depth of box 520 mm  
Width of box 140 mm  
Height of box 146 mm  
Depth of case 539 mm  
Width of case 298 mm  
Height of case 315 mm  
Case volume 0.0506 cbm  
Weight of product 0.0021 kg  
Weight of case 5.4 kg  
EAN of inner box 4066916022914  
EAN case 4066916022921

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2026-03-11

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**