



Product Specification
Replacement filter, for Automatic-Sarpette® , pore
size: 0.2 µm



Product description

Order number	90.189.221
Product description	Replacement filter, for Automatic-Sarpette® , pore size: 0.2 µm, 5 piece(s)/bag

Product characteristics

Suitable for	Automatic-Sarpette®
Pore size	0,2

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	5
Type of smallest subpackaging	bag
Piece(s) / inner box	5
Piece (s) / outer case	500
Piece(s) / pallet	54000

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Replacement filter, for Automatic-Sarpette® , pore
size: 0.2 µm

Depth of case	308 mm
Width of case	223 mm
Height of case	195 mm
Case volume	0.0134 cbm
Weight of product	0.002 kg
Weight of case	1.21 kg
EAN of inner box	4038917419604
EAN case	4038917419598

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

