



Product description

| | |
|---------------------|--|
| Order number | 92.1417.002 |
| Product description | Silicon adapter, for M-Sarpette® , transparent, silicone, 1 piece(s)/bag |

Product characteristics

| | |
|--------------|-------------|
| Suitable for | M-Sarpette® |
|--------------|-------------|

Material & colours

| | |
|-------------------|-------------|
| Product material | silicone |
| Colour of product | transparent |

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Batched | no |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 1 |
| Piece (s) / outer case | 100 |
| Piece(s) / pallet | 32400 |
| Depth of case | 268 mm |
| Width of case | 191 mm |
| Height of case | 100 mm |
| Case volume | 0.0051 cbm |
| Weight of product | 0.006 kg |
| Weight of case | 0.85 kg |
| EAN of inner box | 4038917129008 |
| EAN case | 4038917021906 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

