



Product description

Order number	90.3100.111
Product description	Sarpette® M, single-channel, variable, 1 ml - 10 ml, light blue, 1 piece(s)/case

Product characteristics

Volume of work	1 ml - 10 ml
Coding colour	light blue
Autoclavable (max. °C)	121 °C

Material & colours

Colour of product	grey
-------------------	------

Purity & certification

Product category	no medical device IVD
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Sarpette® M, single-channel, variable, 1 ml - 10 ml, light blue, 1 piece(s)/case

Page 2

Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	1
Depth of box	246 mm
Width of box	88 mm
Height of box	45 mm
Weight of product	0.088 kg
EAN of inner box	4038917937191
EAN case	4038917937207

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2024-03-29

This document was prepared by EDP support and is valid without signature.



SARSTEDT