



### Product description

Order number	90.3100.555
Product description	Sarpette® M, single-channel, variable, 0.5 ml - 5 ml, violet, 1 piece(s)/case

### Product characteristics

Volume of work	0,5 ml - 5 ml
Coding colour	violet
Autoclavable (max. °C)	121 °C

### Material & colours

Colour of product	grey
-------------------	------

### Purity & certification

Product category	no medical device   IVD
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	1
Depth of box	246 mm
Width of box	88 mm
Height of box	45 mm
Weight of product	0.088 kg
EAN of inner box	4038917937146
EAN case	4038917937153

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

