



Product description

Order number	90.3112.200
Product description	Sarpette® M, 12-channel, variable, 20 µl - 200 µl, yellow, 1 piece(s)/case

Product characteristics

Volume of work	20 µl - 200 µl
Coding colour	yellow
Autoclavable (max. °C)	121 °C

Material & colours

Colour of product	grey
-------------------	------

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	1
Depth of box	288 mm
Width of box	150 mm
Height of box	28 mm
Depth of case	290 mm
Width of case	150 mm
Height of case	30 mm
Case volume	0.0013 cbm
Weight of product	0.2 kg
Weight of case	0.52 kg
EAN of inner box	4038917937689
EAN case	4038917937696

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

