




Product Specification
Serological pipette, plugged, 1 ml, sterile, 25
piece(s)/bag

Page 1



Product description

Order number	86.1251.025
Product description	Serological pipette, plugged, max. working volume: 1 ml, coding colour: yellow, material: PS, sterile, 25 piece(s)/bag

Product characteristics

Colour of print/label	yellow
Graduation	yes
Design	plugged
Volume of work	1 ml
Coding colour	yellow

Size

Diameter	4 mm
Length of product	270 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2021-12-05

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
**Serological pipette, plugged, 1 ml, sterile, 25
piece(s)/bag**

Page 2

Material & colours

Product material Polystyrene (PS)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Purity standard sterile
Sterilisation Gamma irradiation
Batched yes

Packaging

Minimum order qty. 1000
Type of inner packaging bag
Piece(s) / inner box 25
Piece(s) / case 1000
Piece(s) / pallet 108000
Case volume 0,0132 cbm
Weight of product 0,0033 kg
Weight of case 3,94 kg
EAN of inner box 4038917127165
EAN case 4038917021074

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2021-12-05

This document was prepared by EDP support and is valid without signature.



SARSTEDT