




Product Specification

**Serological pipette, with tip, plugged, 2 ml, sterile,
non-pyrogenic/endotoxin-free, non-cytotoxic, 1
piece(s)/blister**

Page 1



Product description

Order number	86.1252.001
Product description	Serological pipette, with tip, plugged, 2 ml, coding colour: green, material: PS, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister

Product characteristics

Label/ Print	with printed graduations
Colour of print/label	green/black
Graduation	yes
Design	with tip, plugged
Volume of work	0,2 ml - 2,5 ml
Coding colour	green

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-02-09

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

**Serological pipette, with tip, plugged, 2 ml, sterile,
non-pyrogenic/endotoxin-free, non-cytotoxic, 1
piece(s)/blister**

Page 2

Size

Sample volume	2 ml
Diameter	6 mm
Length of product	270 mm

Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	blister
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	54000
Depth of case	420 mm
Width of case	250 mm
Height of case	290 mm
Case volume	0.0323 cbm
Weight of product	0.0044 kg
Weight of case	5.03 kg
EAN of inner box	4038917127196
EAN case	4038917021081

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-02-09

This document was prepared by EDP support and is valid without signature.



SARSTEDT