




**Product Specification**  
**Serological pipette, with tip, plugged, 2 ml, sterile,  
non-pyrogenic/endotoxin-free, non-cytotoxic, 1  
piece(s)/blister**



**Product description**

Order number 86.1252.001  
Product description Serological pipette, with tip, plugged, 2 ml, coding colour: green, material: PS, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister

**Product characteristics**

Colour of print/label green  
Graduation yes  
Design with tip, plugged  
Volume of work 2 ml  
Coding colour green

**Size**

Diameter 6 mm  
Length of product 270 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

**Serological pipette, with tip, plugged, 2 ml, sterile,  
non-pyrogenic/endotoxin-free, non-cytotoxic, 1  
piece(s)/blister**

**Page 2**

### Material & colours

Product material Polystyrene (PS)  
Colour of product transparent

### Purity & certification

Product category no medical device | IVD  
Purity standard sterile, non-pyrogenic/endotoxin-free, non-cytotoxic  
Sterilisation Gamma irradiation  
Pyrogen-free yes  
Batched yes

### Packaging

Minimum order qty. 1000  
Type of smallest blister  
subpackaging  
Piece(s) / inner box 100  
Piece (s) / outer case 1000  
Piece(s) / pallet 54000  
Depth of case 420 mm  
Width of case 250 mm  
Height of case 290 mm  
Case volume 0.0323 cbm  
Weight of product 0.0044 kg  
Weight of case 5.03 kg  
EAN of inner box 4038917127196  
EAN case 4038917021081

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-06-19

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**